## [EN] Consumer information

**A full contact flexible wound dressing which provides an optimal wound healing environment and abnormal scar prevention.**

Stratamed is a semi-occlusive, self-drying and transparent gel.

When used as directed Stratamed dries to form a protective layer that is gas permeable and waterproof which hydrates and protects compromised skin areas and open wounds from chemical and microbial invasion.

Stratamed can be applied to fresh incisions and excisions, open wounds and compromised skin surfaces including those from post laser treatments, chemical peels, dermal abrasion, 1<sup>st</sup> and 2<sup>nd</sup> degree burns and radiation therapy.

### For open wounds and compromised skin:

Stratamed gel is bacteriostatic and inert.

Stratamed dries to form a thin, flexible wound dressing which hydrates and protects affected areas. Stratamed helps create an optimal wound healing environment for superficial wounds and burns.

This environment leads to faster reepithelialization and a reduced inflammatory response, which reduces recovery time post procedure.

Stratamed does not adhere strongly to the wound bed or the newly formed granulating tissue.

Stratamed is a self-drying, non-sticky, transparent, gel formulation that lightly bonds to the most superficial damaged skin layer.

### For abnormal scar prevention and treatment:

Stratamed softens and flattens raised scars and allows scars to normalize the collagen synthesis process, producing a normal mature scar.

Stratamed can be used with or without a secondary protective dressing.

Stratamed is ideal for exposed areas like the face and neck as well as joints and hairy areas without the need for shaving.

Stratamed is used to relieve itching and discomfort and to reduce redness and discoloration associated with damaged skin and scars.

**Stratamed was developed for use on all types of superficial wounds, burns, skin tears and abrasions including those resulting from:**

- General and cosmetic surgery
- Trauma
- Burns
- Radiation therapy
- Bites
- Grazes
- Laser therapy, including tattoo removal
- Chemical peels
- Dermal abrasion
- Cryotherapy
- Acne and other diseases of the skin

## Intended use

**Stratamed may be applied to open wounds, burns and areas where skin integrity has been compromised.**

### Stratamed is used for:

- The protection and hydration of superficial wounds, burns and compromised skin
- The prevention and treatment of keloids, hypertrophic and excessive scars.

Stratamed may also be used in conjunction with invasive scar management options to improve overall results, including surgical excision and resurfacing, intralesional steroid injections, laser treatments and pressure garments.

Stratamed is suitable for children and people with sensitive skin.

Stratamed may be used on infants, during pregnancy and while breastfeeding.

## Directions for use

Ensure that the affected skin or wound is clean.

Gently pat dry as much excess exudate or wound fluid from the area as possible prior to gel application.

For best results Stratamed should be maintained in continuous contact with the skin (24 hours a day/7 days a week). If Stratamed has been removed by washing it should be reapplied.

When applied correctly to the affected areas, Stratamed should be dry in 5–6 minutes. If it takes longer to dry you have probably applied too much. Gently remove the excess with a clean tissue or gauze and allow the drying process to continue.

### On damaged or broken skin not requiring a secondary dressing:

Apply a very thin layer of Stratamed to the affected area and allow the gel to dry. Stratamed should be applied once daily, or twice daily to exposed areas or as advised by your physician.

Once dry, Stratamed may be covered by sunscreen, cosmetics, pressure garments or casts.

### On damaged skin or wounds requiring a secondary dressing:

Apply a thin layer of Stratamed then cover with the secondary dressing. Drying is not necessary.

Stratamed should be reapplied when changing the dressing or checking wound progress, or as advised by your physician.

## Recommended duration of treatment

Stratamed is recommended for immediate treatment post procedure and should be continued until the skin has fully recovered.

For resurfacing and ablative procedures, continued use is recommended for a minimum treatment of 30 days (24 hours a day/7 days a week) or until no further improvement is seen.

For abnormal scar prevention, Stratamed is recommended

for a minimum treatment of 60–90 days (24 hours a day/7 days a week). Larger scars may take longer.

Continued use is recommended until no further improvement is seen.

## How much Stratamed do I need?

Stratamed gel is a unique formulation that requires substantially less product per application than typical moisturizing creams or barrier ointments.

- **Stratamed 5g (0.17 oz)** is enough to treat a 2–4 cm (1–2 inch) linear incision twice per day for over 90 days.
- **Stratamed 10g (0.35 oz)** is enough to treat a 5–8 cm (2–3 inch) linear incision twice per day for over 90 days.
- **Stratamed 20g (0.7 oz)** is enough to treat a 10–12 cm (4–5 inch) linear incision twice per day for over 90 days or an area of 6x12 cm (2.5x5 inch) twice per day for over 60 days.
- **Stratamed 50g (1.75 oz)** is enough to treat an area of 12x15 cm, (5x6 inch) twice per day for over 60 days.

## Caution

- For external use only.
- Stratamed should not be placed in contact with the eyes.
- Stratamed should not be applied over other skin treatments without the advice of your physician.
- Stratamed may stain clothing if not completely dry.
- Should your wound show signs of infection or fail to heal, consult your physician.
- If irritation occurs, discontinue use and consult your physician.
- Not suitable for highly exudative wounds, tunnelling wounds or 3<sup>rd</sup> degree burns.
- Keep out of the reach of children.
- Store below 25°C (77°F).
- Do not use after the expiration (EXP) date printed on the tube.

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Stratapharma AG