


# Stratamed<sup>®</sup>

advanced film-forming  
wound dressing

## Clearly effective post laser care



The first film-forming silicone gel-based  
wound dressing approved for the use  
on open wounds and compromised skin



CE 1254

Class IIa Medical Device, FDA Listed/TGA approved

 **Stratapharma**  
Switzerland

## A breakthrough in post procedure care

Patients seeking skin rejuvenation and resurfacing treatments have high expectations. They want fast results and the best overall outcome of their therapy. This means reducing the recovery time so patients may return to their daily life with the best possible results.

Following the correct treatment and technique, post procedure care is the key for achieving the best result, patient satisfaction and downtime.

### Use Stratamed immediately post procedure

Stratamed is a full contact flexible wound dressing.

It is the first semi-occlusive, self-drying, transparent silicone wound dressing in a form of gel, that can be applied on compromised and non-epithelialized skin.



- Bacteriostatic, gas permeable and inert<sup>1,2</sup>
- Creates an optimal environment for faster re-epithelialization resulting in reduced downtime<sup>1</sup>
- Reduces and normalizes the inflammatory response<sup>1</sup>
- Hydrates and protects superficial wounds and compromised skin<sup>1</sup>
- Does not adhere strongly to the wound bed or the newly formed granulating tissue<sup>2</sup>
- Improves the visible outcome of the treatment<sup>1</sup>
- Provides symptomatic relief (redness/discoloration, itching, discomfort, pain)<sup>1</sup>
- Easy to apply with a great safety profile<sup>1</sup>



Day 0



Day 7

#### Laser resurfacing (acne scars)<sup>1</sup>

During treatment and 7 days post treatment with Stratamed in combination with thermage and CO<sub>2</sub> fractional laser resurfacing.

*"We observed a faster re-epithelialization of the wound surface with visible satisfactory results starting as early as 7 days, compared with that we observed with our standard practice"<sup>1</sup>*

**Stratamed can be used immediately after all dermatological interventions like:**  
Laser treatment, dermal abrasion, chemical peels, tattoo removal, surgery.

# Clinically proven to be effective and easy to use<sup>1-3</sup>



**Laser resurfacing (wrinkles)<sup>1</sup>**  
 Before, during and 10 days post treatment with Stratamed in combination with thermage, CO<sub>2</sub> fractional laser resurfacing and Botox.

Before treatment      Day 0      Day 10



**Dermabrasion<sup>1</sup>**  
 During treatment and 4 weeks post treatment with Stratamed (monotherapy).

During treatment      Day 0      Week 4



**Rapid healing of chronic wound resulting from surgery after squamous cell carcinoma<sup>2</sup>**  
 Nonhealing wound before (4 months post surgery) and 7 days post treatment with Stratamed (monotherapy).

4 months post surgery      Day 7

Stratamed **speeds up** the healing process resulting in reduced downtime.<sup>3</sup>  
 Stratamed **improves** the visible outcome of the treatment.<sup>3</sup>

Stratamed in monotherapy		Thermal water + white petrolatum	
			
Immediately after ablative fractional resurfacing	7 days post treatment with Stratamed	Immediately post procedure	15 min after primary dressing application
			
Immediately after ablative fractional resurfacing	7 days post treatment with thermal water + white petrolatum		

Acute erythema measured with hemoglobin variation 15 min after application

<b>+1%</b>	Thermal water + white petrolatum
<b>-31%</b>	Stratamed monotherapy

15 min after Stratamed application

Antera 3D measurement

Stratamed dries into a full-contact primary dressing over microtextural surface irregularities - excellent soothing effect.<sup>3</sup>

Stratamed immediately decreases post-inflammatory burning sensation by reducing erythema and superficial skin temperature.<sup>3</sup>

Stratamed leads to a more comfortable post procedure care and increases user's adherence to multiple session protocols.<sup>3</sup>



Stratamed is available in 5g, 10g, 20g and 50g tubes.

# Stratamed®


advanced film-forming wound dressing

## Directions for use

For best results Stratamed should be maintained in continuous contact with the skin (24 hours a day / 7 days a week). When applied correctly to exposed areas, Stratamed should be dry in 5–6 minutes. If it takes longer to dry you have probably applied too much. Gently remove the excess with a clean tissue or gauze and allow the drying process to continue.


- Ensure that the affected skin or wound is clean.
- Pat dry as much excess exudate or wound fluid from the wound as possible prior to gel application.

### Post skin rejuvenation and resurfacing treatments:



**On damaged or broken skin not requiring a protective secondary dressing:**  
Apply a very thin layer of Stratamed to the affected area and allow the gel to dry. Stratamed should be applied once daily, or twice daily to exposed areas or as advised by your physician. Once dry, Stratamed can be covered by sunscreen, cosmetics, pressure garments or casts.

### After surgical incisions/excisions:



**On damaged skin or wounds requiring a protective secondary dressing:**  
Apply a thin layer of Stratamed then cover with the secondary dressing. Drying is not necessary. Stratamed should be reapplied when changing the dressing or checking wound progress, or as advised by your physician.

## From Stratpharma Switzerland, producer of:

### Stratamark®

for the prevention and treatment of stretch marks

**The novel formulation for the prevention and treatment of Striae.**

Can be used in combination with non ablative lasers for improving overall results.



[www.stratamark.net](http://www.stratamark.net)



### Strataderm®

scar therapy gel

**For the prevention and treatment of hypertrophic scars and keloids.**

For treatment of old and new scars. Softens and flattens raised scars.



[www.strataderm.com](http://www.strataderm.com)



[www.stratamed.com](http://www.stratamed.com)

References: 1. Sandhofer M, Schauer P. SKINmed 2012;10:51–57  
2. Monk EC, Benedetto EA, Benedetto AV. Dermatologic Surgery 2014;40:76–79  
3. Data on file. Stratpharma AG

